

***IN THE UNITED STATES PATENT AND TRADEMARK OFFICE***

Applicant: H. William Bosch et al.  
Title: NOVEL NIMESULIDE COMPOSITIONS  
Appl. No.: 10/697,703  
Filing Date: 10/31/2003  
Examiner: Tristan J. MAHYERA  
Art Unit: 1615  
Confirmation Number: 8369

DECLARATION UNDER 37 CFR 1.131

Commissioner for Patents  
P.O. Box 1450  
Alexandria, VA 22313-1450

Sir:

I, H. William Bosch, hereby declare and state that:

1. I am a citizen of the United States, residing at 237 Rodney Circle, Bryn Mawr, PA 19010.
2. At the time of events detailed in paragraph 4, *infra*, I was an employee of Elan Drug Delivery, Inc., with offices at 3500 Horizon Drive, King of Prussia, PA 19406.
3. I am a co-inventor of the invention disclosed and claimed in the above-referenced application.

4. Prior to June 27, 2003, I instructed my associates, as part of my supervisory role, to prepare nimesulide compositions comprising particles of nimesulide or a salt thereof having an effective average particle size of less than 2000 nm and at least one surface stabilizer adsorbed on the surface of the particles. My work relating to preparing the nimesulide compositions, which occurred prior to June 27, 2003, is documented in the attached exhibits.

5. As shown in Exhibit A (Notebook No. 5822, pages 006-008), the formulation comprising 5% nimesulide and 1% Plasdone<sup>®</sup> S-630 provides a stable nanoparticulate nimesulide composition.

6. As shown in Exhibit B (Notebook No. 5822, pages 009-011), the formulation comprising 5% nimesulide and 1% Plasdone<sup>®</sup> S-630 provides a stable nanoparticulate nimesulide composition.

7. As shown in Exhibit C (Notebook No. 5822, pages 012-014), the formulation comprising 5% nimesulide, 1% Plasdone<sup>®</sup> S-630 and 0.2% DOSS provides a stable nanoparticulate nimesulide composition.

8. As shown in Exhibit D (Notebook No. 5822, pages 015-017), the formulation comprising 5% nimesulide, 1% Plasdone<sup>®</sup> S-630 and 0.05% sodium lauryl sulfate (SLS) provides a stable nanoparticulate nimesulide composition.

9. I further declare that all statements made herein of my own knowledge are true and that all statements made on information and belief are believed to be true; and further, that these statements were made with the knowledge that willful false statements and the like so made are punishable by fine or imprisonment, or both, under Section 1001 of Title 18 of the United States Code and that such willful false statements may jeopardize the validity of the application or any patent resulting therefrom.

Date

June 27, 2008

H. William Bosch

H. William Bosch

Title Nimesulide (5% API ; 1% S-630)(cont. from pg. 005)**Batch Record for Dispersion Technology Milling Procedures****I. General Information**

Name	<u>Christian Wertz</u>
Date	
Formula	<u>5% - Nimesulide 1% S-630</u>
Continued on Page	<u>007</u>

**II. Quantities Dispensed**

	Quantity	Type	Source	Lot Number
Media	<u>80.6</u>	<u>Polymill 200</u>	<u>Dow / PMRS</u>	
Drug Substance	<u>4.25</u>	<u>Nimesulide</u>	<u>Sigma</u>	
Stabilizer	<u>0.85</u>	<u>S-630</u>		
Water	<u>79.9</u>	<u>DI</u>		
Other				

**III. Process Parameters**

Milling Method	<u>Dynomill (150 cc batch chamber) F915 @ Rm. 205</u>
Mill Speed	<u>4,200 rpm</u>
Temperature	<u>~10°C</u>

**IV. Notes**

Milling Time:	<u>9:53 - Start milling ; 10:39 - first sample</u>
	<u>10:53 - second sample ; 11:41 - third sample</u>
	<u>1:10 - Harvest</u>
Quantity retained post-milling:	<u>forget to filter out media before weighing</u>

\* Did not filter out media initially and discarded ~ 1/2 suspension.  
Later filtration left ~20 mL of media free suspension.

**CONFIDENTIAL**

Signature

Christian F. Wertz

Date

Reviewed and understood by

Klaus Riedel

Date

(cont. on pg. 007)



Date \_\_\_\_\_





Title Nimesulide (5 % API, 2 % S-630)

(cont. from pg. \_\_\_\_\_)

- mix S-630 slowly into DI H<sub>2</sub>O w/ mild stirring until dissolved
- add polymill 200 w/ gentle manual stirring
- add API w/ gentle stirring until thoroughly mixed

### Batch Record for Dispersion Technology Milling Procedures

#### I. General Information

Name	Christian Wertz
Date	
Formula	5 % Nimesulide, 2 % S-630
Continued on Page	010

#### II. Quantities Dispensed

	Quantity	Type	Source	Lot Number
Media	80.6	Polymill 200	DOW	MM001012
Drug Substance	4.25	Nimesulide	Sigma	117H1019
Stabilizer	1.70	S-630	ISP Tech.	ML90012974
Water	79.05	DI		
Other				

#### III. Process Parameters

Milling Method	Dynomill (150 cc batch chamber) F915 @ Rm. 205
Mill Speed	4200 rpm
Temperature	~10 °C

#### IV. Notes

Milling Time:	8:15 Start batch ; 9:15 1 <sup>st</sup> Sample
	10:15 2 <sup>nd</sup> Sample ; 11:15 Harvest
Quantity retained post-milling:	49.6 g (58 %)

(cont. on pg. 010)

**CONFIDENTIAL**

Signature

Christian F. Wertz

Date

Reviewed and understood by

H. William Bosch

Date







EXHIBIT C



Title Nimesulide 5% API, 1% S-630, 0.2% DOSS(cont. from pg. 013 )**Particle Size Stability for Dispersion Technology Formulations****I. General Information**

Name	Christian F. Wertz
Date	
Formulation	5% API, 1% S-630, 0.2% DOSS
Notebook reference	
Continued on page	

**II. Particle Size Data**

Particle Size Analyzer Used	HORIBA LA-910 (s#: 8514870103D)
Standards Measured	Lot #: 22569; mean = 203 ; Duke Sci.; 200 nm standard

Elapsed Time	Storage Conditions	Mean, nm	D50, nm	D90, nm	Comments
1 day	5°C	136	116	223	no Sonication
	5°C	137	116	224	60 s Sonication
2 day	5°C	143	121	238	no Sonication
	5°C	144	122	241	60 s Sonication
5 day	5°C	149	133	239	no Sonication
	5°C	151	135	242	60 s Sonication
7 day	5°C	160	143	259	no Sonication
	5°C	163	146	261	60 s Sonication
21 day	5°C	162	150	252	no Sonication
	5°C	166	155	255	60 s Sonication
35 day	5°C	180	172	276	no Sonication
	5°C	187	180	280	60 s Sonication

<sup>cfw</sup>  
Data in folder ~~to~~ Supplementary folder CFW-5822 A

**CONFIDENTIAL**

(cont. on pg. \_\_\_\_\_)

Signature

Christian F. Wertz

Date

Reviewed and understood by

Ad. Wilhelm Bosch

Date



Title Nimesulide 5% API, 1% S-630, 0.05% SLS

(cont. from pg. \_\_\_\_\_)

- Dissolved S-630 in H<sub>2</sub>O followed by SLS under gentle mixing
- SLS dissolved very rapidly w/ very little foam

### Batch Record for Dispersion Technology Milling Procedures

#### I. General Information

Name	Christian F. Wertz
Date	
Formula	5% API, 1% S-630, 0.05% SLS
Continued on Page	

#### 4. Quantities Dispensed

	Quantity	Type	Source	Lot Number
Media	80.6	Polymill 200	DOW	MM001012
Drug Substance	4.25	Nimesulide	Sigma	117H1019
Stabilizer	0.85	S-630	ISP Technology	ML900012974
Water	79.86	H <sub>2</sub> O	DI	
Other	0.04	SLS		

#### III. Process Parameters

Milling Method	Dynomill (150 cc batch chamber) F915 @ room 205
Mill Speed	4200 rpm
Temperature	10 C

#### IV. Notes

Milling Time: 8:22 Start batch ; 9:22 1 <sup>st</sup> sample
10:22 Harvest
Quantity retained post-milling: 69.4 g (80.7 wt %)

- mill began leaking after first sample was taken from mill

(cont. on pg. 016)

**CONFIDENTIAL**

Signature

Christian F. Wertz

Date

Reviewed and understood by

W. Wilhelm Bosch

Date



Title Nimesulide 5% API, 1% S-630, 0.05% SLS(cont. from pg. 016)**Particle Size Stability for Dispersion Technology Formulations****I. General Information**

Name	Christian F. Wertz
Date	
Formulation	5% API, 1% S-630, 0.05% SLS
Notebook reference	
Continued on page	

**II. Particle Size Data**

Particle Size Analyzer Used	HORIBA LA-910 (s#: 8514870103D)
Standards Measured	Lot #: 22569; mean = ; Duke Sci.; 200 nm standard

Elapsed Time	Storage Conditions	Mean, nm	D50, nm	D90, nm	Comments
3 day	5°C	123	108	192	no sonication
	5°C	123	108	193	60 s sonication
5 day	5°C	127	110	203	no sonication
	5°C	127	110	205	60 s sonication
7 day	5°C	133	113	219	no sonication
	5°C	134	113	219	60 s sonication
21 day	5°C	126	110	199	no sonication
	5°C	127	109	205	60 s sonication
35 day	5°C	141	120	234	no sonication
	5°C	142	119	237	60 s sonication

Data in folder #2 supplementary folder CFW-5822 A

**CONFIDENTIAL**

Signature

Christian F. Wertz

Date

Reviewed and understood by

Ad. William Bosch

Date

(cont. on pg. \_\_\_\_\_)